

Excellent

ENT treatment unit



Ergonomic access
to all your therapeutic equipment

MediTop | Excellent

The Excellent unit offers a design that optimally uses its space and is ergonomic and hygienic. This unit differs itself from other models due to the ease of installation and the extremely wide variety of customization and is therefore very suited for all customers.

Optimal

This highly requested and totally equipped unit comes in 3 models. As a result, it adapts perfectly to everyone's wishes. Almost all additional functions can be fully integrated in the Excellent units, allowing the

It adapts perfectly to everyone's wishes

doctor to work very smoothly and ergonomically. The unit also offers a lot of storage space for instruments and disposable material.


Hygienic

Due to the use of easy-to-maintain materials of high quality and its easy-to-clean sleek design, these units score very high in terms of hygiene. The units are virtually seamless and the absence of wall and floor gaps allows for a quick and thorough cleaning of the unit with the most commonly used household products.

Customization

The Excellent models are designed in such a way that the provided spaces can be used in an optimal way.





The units can all be expanded with more blocks, that add new functions to the base unit. Usually a sanitary block is added. With the extra additions, the worktops also expand and thus provide more work space.

Wide range of material finishes

This model is available with a Marble Top finish or with HPL (High Pressure Laminate). When choosing HPL, a whole range of design possibilities such as wood/metal/stone looks become available.

The provided spaces can be used in an optimal way

A wide selection of integrated functions

The ExcellENT can be equipped with a wide selection of integrated functions to accommodate specific requirements.

Standard integrated functions

- 2 instrument trays and 4/5 handpiece holders
- 1 large storage-drawer
- 1 drawer for hygienic storage of used instruments and waste

Integrated functions by choice

- Suction system (variants available)
- Compressed air system (variants available)
- LED light source (variants available)
- Ear syringe (AquaTop XS)
- Cautery system (complete with handpiece and electrodes)
- Mirror heater (pre-heater and/or hot air)

- Rigid Scope holders (heated or non-heated)
- Lateral spittoon system
- Video-endoscopy
- Video-stroboscopy
- Microscope with integrated camera
- HF-surgery / coagulation device
- Tympanometer
- Rhinomanometer
- Echo / sinus A-scan module
- Image & Video Manager System (IVM)
- Hygiene tracker
- Integration of sanitary blocks, extra cupboards, custom-made desktops in the same material and colour as your treatment unit, etc...
- Chair control
- Light control (DALI, etc)
- Blinds control

MediTop | ExcellENT

The ExcellENT family

ExcellENT Microscopy

The ExcellENT Microscopy units are specially designed for treatments in mainly otologic settings, where patients are examined while lying down.

To facilitate the examinations, the microscopy units are equipped with 2 instrument parts; 1 on each side of the extended therapy section. These units can be equipped with many more options, which makes working with these units even more pleasant.



ExcellENT II

The ExcellENT II unit is a simplified version of the ExcellENT unit, ideal for settings where a high degree of customization and integration of options is not required.

This unit nevertheless offers the possibility to integrate the most requested options and thus allows smooth and ergonomic work. The ExcellENT II unit has a smaller touchscreen, a thinner top and only 4 handpiece holders instead of 5.

ExcellENT IIS





The ExcellENT IIS unit is the smallest unit within the ExcellENT group, making the integration of extra options more limited.

Despite this limitation, this unit also allows smooth and ergonomic work.



MediTop | Excellent


Standard dimensions



Contact us
sales@meditop.com

We would like to meet you and your company
and see how we can help you
optimize your work environment.

Optional accessories



AQUATOP XS | Ear syringe



MEDITOP LED HEADLIGHT



MEDITOP COMFORT CHAIR




IMAGE & VIDEO SOFTWARE



ZEISS | microscopes

Technical specifications

Voltage	230V
Power	2300VA
Dimensions	3150 x 600 x 950 mm (incl sanitary part)