

MT750

Mobile ENT treatment unit



Our smallest mobile ENT treatment unit

Product code: 2000350/2000500/2000550
MT750

MediTop I MT750

The MT750 ENT unit is the smallest ENT unit MediTop produces. This elegant model is equipped with 3 instrument holders for suction, compressed air, cold light system etc. It is a perfect mobile solution for small practices and the A&E department.

Mobility

The MT750 is highly moveable due to its 4 castors. Despite the unit being easy to move around, it is robust and reliable. The unit only needs one electrical power socket to be able to function.

Highly moveable ENT treatment unit

Simplicity


The MT750 ENT unit accommodates easy operating of its integrated functions. There are integrated functions

that start and stop automatically by taking out or placing back the handpieces in the handpiece holders. Other functions can be easily activated by the push of a button.

Customizable

Every MT750 ENT unit is built to order. That makes it possible to produce the unit in accordance with customer requirements when it comes to materials and finish. This allows for specific colour coordination in the consulting room. There is an almost endless range of colours available.





- A Apparatus compartment
- B Handpiece holders
- C Instruments compartment
- D Waste bin
- E Used instruments
- F Electronics compartment
- G Storage for consumables
- H Main power switch (side of the unit)
- J Control panel

The MT750 in more detail

- A. Several small devices are placed in this compartment.
- B. The corresponding integrated functions are activated automatically whenever the handpieces are removed from their holders.
- C. This compartment houses several instrument trays.
- D/E. The waste (D) and the used instruments (E) can be placed in two separate containers.
- F. Behind this panel you can find the electronics compartment. Also located here is the on/off switch.
- G. The drawer (G) offers space for instruments and consumables.
- H. With the main power switch, the complete MT750 basic treatment system can be controlled (on or off).
- J. Main control panel.

Basic unit with all necessary functions

A wide selection of integrated functions

The MT750 can be equipped with a wide selection of integrated functions to accommodate specific requirements.

Standard integrated functions

- 2 instrument trays and 3 handpiece holders
- 1 large drawer for hygienic storage of used instruments and waste
- 1 large storage-drawer


Integrated functions by choice

- Suction system (variants available)
- Compressed air system (variants available)
- LED light source (variants available)
- Mirror heater (pre-heater and/or hot air)
- HF-surgery / coagulation device

Contact us
sales@meditop.com

We would like to meet you and your company to see how we can help you optimize your work environment.


Optional accessories



AQUATOP XS | Ear syringe



MEDITOP LED HEADLIGHT



MEDITOP COMFORT CHAIR



MICROSCOPES



ENDOSCOPES

Technical specifications

Voltage	230V
Power	2300VA
Dimensions (L x D x H)	991 x 575 x 935 mm