

MediTop Comfort

Safe, ergonomic & comfortable patient chair
with extra low access




The MediTop Comfort chair distinguishes itself due to the high level of comfort for both patient and doctor and the integration with MediTop ENT treatment units.

It is designed with an extra low level of access making it easy for all patients to sit in it. The use of high quality motors and a smart design ensure quiet and smooth operation.



Benefits at a glance

- Electric height adjustment between 49-79 cm
- Multi-functional footswitch
- Automatic return
- Electric backrest reclines between 90° and 180°
- Adjustable headrest
- Rotation up to 340°
- Footrest can be folded away
- Electric leg support moves in sync with the backrest
- Many colours available
- Options: mobile version available, paper roll holder available, automatic height adjustment by facial recognition, extended leg support, electric brake function, integrated wheels
- Also available: MediTop Basic: backrest with gas spring control



Technical specifications

Width seat	49 cm
Seat height (min/max)	49/79 cm
Total width	71 cm (including armrest)
Length lying surface	182 cm
Chair	340° rotatable
Maximum power use	120 VA
Voltage	230 Volt AC - 50/60 Hz
Maximum weight	150 kg