

## MediTop LED headlight

A lightweight, cordless headlight for outpatients and theatre use



The MediTop LED headlight is a comfortable lightweight, cordless examination headlight that is suitable for prolonged use in the OPD and OT. The headlight automatically switches on whenever it is removed from its holder and it switches off when placed back on the holder. This saves the user from having to switch the headlight on and off.

Three light intensity settings are controlled by the button on the headlight.



## Benefits at a glance

- Light weight
- Comfortable, adjustable headband
- 3 light intensity settings
- Automatic off switch (motion sensor)
- Stand-alone or wall holder
- Automatic charging with supplied holder
- Battery charger and holder can be integrated in ENT units
- Simple intuitive control
- Remaining charge LED
- Long battery life



## Technical specifications

Power adapter	12V – 2.5A
Weight	290 gram
Battery life	4.5 hours
Colour temperature	5000K
Diameter light spot at 420 mm	120 mm