

Innovia Medical®

Nasal Splints & Packing



**SAMPLES
AVAILABLE**

and NEW 2024
price list

INTERNAL & EXTERNAL

KEY PRODUCT FEATURES

THE DENVER SPLINT®


- ✓ Moldable aluminum
- ✓ No water mess
- ✓ Water Resistant

BREEZE™ NASAL SPLINTS

- ✓ Unique airway geometry creates a pressure gradient force that increases airflow
- ✓ Suture indicators are located at anterior tip to aid suture placement. Surface characteristics aid insertion and removal.

RHINO ROCKET®

- ✓ Sterile PVA foam
- ✓ Epistaxis management after trauma, nasal or sinus surgery
- ✓ Gentle but firm mucosal compression



CONTENTS

INTERNAL & EXTERNAL NASAL SPLINTS & NASAL PACKING

THE DENVER SPLINT® → Page 4

THE DENVER SPLINT® THERMOPLASTIC NASAL SPLINTS → Page 6

INTERNAL NASAL SPLINTS → Page 7

NASAL SEPTAL BUTTONS
BREEZE™ NASAL SPLINTS
BIVALVE NASAL SPLINTS

NASAL AND SINUS PACKS → Page 8

ABSORBENT™ NASAL AND SINUS PACKS

RHINO ROCKET® → Page 9

EPISTAX® EPISTAXIS BALLOON CATHETERS → Page 9

SAMPLES AVAILABLE ON REQUEST

Distributed by DTR Medical Ltd (an Innovia Medical Company)

17 Clarion Court, Clarion Court Enterprise Park, Swansea, SA6 8RF, UK.


Please work with your Business Manager for samples, or contact us in one of the following ways:

Call: +44 (0) 1792 79 79 10 **Email:** dtrinfo@innoviamedical.com

THE DENVER SPLINT®

THE DENVER SPLINT®


THE DENVER SPLINT® - TWO PIECE



The Denver Splint® Series 1500

- Soft moldable aluminum and Velcro™
- Length of wear 7-10 days

PART NO.	DESCRIPTION	DIMENSIONS	QTY/BOX
10-1500-10NKM	Mini	50mm x 28mm	10/box
10-1500-10NKP	Petite	50mm x 30mm	10/box
10-1500-10KNS	Small/Medium	65mm x 35mm	10/box
10-1500-10NKL	Large	69mm x 43mm	10/box



The Denver Splint® Series 5000

- Soft moldable aluminum and Velcro™ with attached dorsal pad
- Length of wear 7-10 days

PART NO.	DESCRIPTION	DIMENSIONS	QTY/BOX
10-5000-10NKM	Mini	50mm x 28mm	10
10-5000-10NKP	Petite	50mm x 30mm	10
10-5000-10NKS	Small/Medium	65mm x 35mm	10
10-5000-10NKL	Large	69mm x 43mm	10

THE DENVER SPLINT® - ONE PIECE



The Denver Splint® Series 2000


- Soft moldable aluminum and 1/8" ventilating foam
- Length of wear 5-7 days

PART NO.	DESCRIPTION	DIMENSIONS	QTY/BOX
10-2000-10NKM	Mini	50mm x 28mm	10
10-2000-10NKP	Petite	50mm x 30mm	10
10-2000-10NKS	Small/Medium	65mm x 35mm	10
10-2000-10NKL	Large	69mm x 43mm	10

SPLINTING STEPS



1. Dorsal Pad



2. Adhesive Velcro



3. Aluminum Splint



- Adhesive
- Moldable aluminum
- No water mess
- Inconspicuous
- Water resistant

- Stabilization of the nose after surgery or trauma
- Provides compression to aid the reduction of post-operative edema
- Moldable aluminum retains nasal bones in position
- Layers give cushioned pressure on subcutaneous tissue
- Dorsal pad reduces dead space and possibility of pollybeak deformity
- Four sizes: Mini, Petite, Small/Medium, Large

Adhesive | Moldable aluminum | No water mess

FREE SAMPLES ON ALL PRODUCT CODES!

THE DENVER SPLINT® THERMOPLASTIC NASAL SPLINTS

THE DENVER SPLINT® THERMOPLASTIC NASAL SPLINT

- Rigid, opaque, lightweight material
- Soften with hot water and mold over taped nose

Custom Cut Options Available:


- 3x3" Squares
- 9x12" Sheets



THE DENVER SPLINT® - SERIES 4000 (PRE-CUT)

PART NO.	SIZE	DESCRIPTION	QTY/BOX
10-4021-10NKL	Large	Ivory 1/16"	10
10-4021-10NKM	Mini	Ivory 1/16"	10
10-4021-10NKP	Petite	Ivory 1/16"	10
10-4021-10NKS	Small	Ivory 1/16"	10
10-4022-10NKL	Large	Blush 1/16"	10
10-4022-10NKM	Mini	Blush 1/16"	10
10-4022-10NKP	Petite	Blush 1/16"	10
10-4022-10NKS	Small	Blush 1/16"	10

Indicate Size:



INTERNAL NASAL SPLINTS

BREEZE™ NASAL SPLINTS

Unique airway geometry creates a pressure gradient force that increases airflow. Suture indicators are located at anterior tip to aid suture placement. Surface characteristics aid insertion and removal. Sterile, single-use.

- Increase airflow
- Soft medical-grade silicone
- Commonly referred to as "Doyle Splint"

PART NO.	DESCRIPTION	QTY
SP-79300	With Airway, Sterile	1 Pair
SP-79300-05	With Airway, Sterile	5 Pairs
SP-79300-50	With Airway, Sterile	50 Pairs
SP-79400	Without Airway, Sterile	1 Pair



BIVALVE NASAL SPLINTS




Non-stick surface reduces adhesion and/or hematoma. Medial slit facilitates insertion and removal. Holes for easier suturing. Sterile, single use. Clear fluoroplastic, medical-grade material.

PART NO.	DESCRIPTION	QTY
SP-27000	Standard Thin (.25mm)	6 Pairs/box
SP-27005	Standard Thick (.5mm)	6 Pairs/box
SP-27010	Oversized Thin (.25mm)	6 Pairs/box
SP-27015	Oversized Thick (.5mm)	6 Pairs/box

FREE SAMPLES ON ALL PRODUCT CODES -
CONTACT YOUR BUSINESS MANAGER TODAY!

NASAL AND SINUS PACKS


ABSORBENT™ NASAL AND SINUS PACKS



AbsorbENT™ nasal and sinus packs utilize three-dimensional interconnected open cells that are soft and pliable.


- Even and fast expansion
- Soft, lint and fiber free
- Rapid wicking properties

Note : Dimensions listed are for hydrated state.



PART NO.	SIZE	DIMENSIONS LxWxH (CM)	QTY
RH-2315-20	Small	5.5 x 1.5 x 2	20
RH-7401-10	Small	4.5 x 1.5 x 2	10
RH-7408-10	Small with Airway	4.5 x 1.5 x 2	10
RH-7404-10	Large	8 x 1.5 x 2	10
RH-7405-10	Large with Airway	8 x 1.5 x 2	10
RH-7409-10	Large Anatomical	8 x 1.5 x 3	10
RH-7410-10	Large Anatomical with Airway	8 x 1.5 x 3	10
RH-7411-10	Large Anatomical - Thin	8 x 1 x 3	10
RH-7406-10	Large - Thin	8 x 1 x 2	10
RH-2310-10	Extra Large	10 x 1.5 x 2.5	10

SINUS PACKS



PART NO.	DESCRIPTION	DIMENSIONS LxWxH (CM)	QTY
RH-7530-10	Sinus, 6W	3.5 x 1.2 x 0.6	10
RH-7532-10	Sinus, 6W	3.5 x 1.2 x 0.9	10

SAMPLES AVAILABLE

New 2024 price list

Please work with your Business Manager for samples, or contact us in one of the following ways:

Call: +44 (0) 1792 79 79 10 **Email:** dtrinfo@innoviamedical.com



RHINO ROCKET®

- Sterile PVA foam
- Epistaxis management after trauma, nasal or sinus surgery
- Gentle but firm mucosal compression
- Tamponade of anterior or middle segment epistaxis
- 30 Seconds to position and insert
- Foam swells to six times its compressed diameter
- String attached for retention and ease of removal
- Dimensions listed are for hydrated state and for insertion position

SLIMLINE



PART NO.	SIZE	DIMENSIONS LxWxH (CM)	QTY
11S-S0300-08AS	Small	3 x 1 x 2	8
11S-S0500-08AS	Medium	5 x 1 x 2	8
11S-S0800-08AS	Large	8 x 1 x 2	8
11S-S1000-08AS	Extra Large	10 x 1 x 2	8

WIDELINE

PART NO.	SIZE	DIMENSIONS LxWxH (CM)	QTY
11W-W0300-08AS	Small	3 x 1.5 x 2	8
11W-W0500-08AS	Medium	5 x 1.5 x 2	8
11W-W0800-08AS	Large	8 x 1.5 x 2	8
11W-W1000-08AS	Extra Large	10 x 1.5 x 2	8

STARTER-PACK

PART NO.	DESCRIPTION	QTY
11S-STARTERRHINO	Variety pack with one of each size.	8/box



STARTER-PACK

PART NO.	DESCRIPTION	QTY
11S-STARTERRHINO	Variety pack with one of each size.	8/box

EPISTAX® EPISTAXIS BALLOON CATHETERS


Double balloon catheter made from soft, medical-grade silicone, designed to control intranasal bleeding. Allows for control of anterior and/or posterior nasal bleeding. Available with or without Syringe. Made from soft, medical-grade silicone. Sterile. Single-use.



PART NO.	SIZE	DESCRIPTION	QTY/BOX
RH-4200-01	97mm	No Syringe	1
RH-4200-10	97mm	No Syringe	10
RH-4220-01	97mm	With 20cc Syringe	1
RH-4400-01	102mm	No Syringe	1
RH-4400-10	102mm	No Syringe	10
RH-4420-01	102mm	With 20cc Syringe	1

INNOVIA MEDICAL

We also
manufacture
a range of
single-use
ENT devices
that compliment
our Rhinology range



RESPONSIBLE MANUFACTURERS OF MEDICAL DEVICES

As part of Innovia Medical's effort to be a more sustainable and environmentally friendly company, we have begun our Manufacturing journey to becoming net zero by 2040. We believe in protecting the planet by reducing our carbon footprint.

For samples contact us in one of the following ways:

Call: +44 (0) 1792 79 79 10

Email: dtrinfo@innoviamedical.com

www.innoviamedical.com

innoviamedical.com/net-zero



Innovia Medical

Nasal Splints & Packing

Manufactured by

Summit Medical, LLC (an Innovia Medical Company)

815 Vikings Parkway, Suite 100. St. Paul, MN 55121

1-888-229-2875

www.innoviamedical.com

Distributed by

DTR Medical Ltd. (an Innovia Medical Company)

17 Clarion Court, Clarion Court Enterprise Park, Swansea, SA6 8RF, UK.

Call: +44 (0) 1792 79 79 10

Email: dtrinfo@innoviamedical.com